

# **CLINICAL APPLPLICATIONS OF NEXT GENERATION SEQUENCING**

Bharat Thyagarajan

Department of Laboratory Medicine and Pathology

University of Minnesota

# MOLECULAR DIAGNOSTICS LABORATORY

- The Molecular Diagnostics Laboratory (MDL) processes around 20,000 specimens annually
- Major testing categories
  - Infectious disease testing: HPV
  - Bone marrow engraftment analyses
  - Hematological malignancies
    - Translocations, quantitative BCR-ABL, JAK2/FLT3/NPM1/CEBPA. T and B cell gene rearrangements etc.
  - Solid tumor malignancies
    - Microsatellite instability, KRAS, BRAF etc.
  - Inherited disorders
    - Factor II, V mutations, sequencing, Southern blot etc.

# CURRENT MOLECULAR TESTING SCHEME FOR ONCOLOGY

CLINICIANS ORDER INDIVIDUAL GENETIC TESTS

MOLECULAR DIAGNOSTICS

DNA/RNA EXTRACTED FROM  
SUBMITTED TISSUE

ONE MUTATION = ONE TEST

SEPARATE MOLECULAR PATHOLOGY REPORT

CYTOGENETICS

FISH

KARYOTYPE

ARRAY CGH

SEPARATE CYTOGENETICS REPORTS


# LUNG CANCER

- As recently as a decade ago stage IV lung cancer had an universally poor prognosis of < 12 months irrespective of chemotherapeutic regimen
- At present, lung cancer with EGFR mutations have a mean survival of > 2 years
- Initial results of targeted therapy for other genetic alterations (e.g. ALK, ROS translocations) have shown promise
- Thus lung cancer is now considered a prototype for genetically tailored cancer therapy

# CURRENT GENETIC TESTING FOR LUNG CANCER

- Molecular Diagnostics
  - *EGFR* mutation analysis
    - Wide range of genetic alterations including point mutations in various exons (18-21) and deletions in exon 19
- Cytogenetics
  - *ALK-EML4* translocations
    - Commonly detected using an *ALK* break-apart probe

# MUTATIONAL PROFILE IN LUNG CANCER



# CURRENT MOLECULAR TESTING SCHEME FOR INHERITED DISEASE

CLINICIANS ORDER INDIVIDUAL GENETIC TESTS

MOLECULAR DIAGNOSTICS

DNA EXTRACTED FROM BLOOD

ONE GENE = ONE TEST

EACH TEST IS A SEPARATE MOLECULAR  
PATHOLOGY REPORT


CYTOGENETICS

FISH


KARYOTYPE

ARRAY CGH

SEPARATE CYTOGENETICS REPORTS




# DISTRIBUTION OF MSI vs. MSS COLON CANCERS






# TESTING FOR LYNCH SYNDROME


## Familial Influences in Colorectal Cancer



## Colorectal Carcinoma Relative Frequency by Pathogenesis



# CURRENT TESTING ALGORITHM



# **LIMITATIONS OF CURRENT TESTING PARADIGM**


- **INHERITED DISORDERS**
  - Comprehensive genetic testing for several syndromes frequently involve simultaneous testing for several genes
  - Increasing demand for detection of point mutations and structural genetic alterations within tested genes
- **CANCER DIAGNOSTICS**
  - Comprehensive prognostic and predictive testing in near future will involve testing at least a few dozen genes
  - Various types of genetic alterations (point mutations, translocations etc.) will need to be evaluated simultaneously
  - Limited amount of sample available will be available for testing

# PROPOSED SOLUTION


- NEXT GENERATION SEQUENCING technology was specifically designed to simultaneously evaluate variation in several genes
- This technology can also be used to detect different types of genetic alterations
- TYPES OF SEQUENCERS
  - HiSeq 2000/2500
  - Desktop sequencers: MiSeq/IonTorrent



# TECHNOLOGY: NEXT GENERATION SEQUENCING



# WHOLE GENOME VS. TARGETED CAPTURE



# MAJOR STEPS OF NGS


- DNA library preparation
- Target enrichment
- Cluster generation & sequencing (Illumina HiSeq 2000)
- Bioinformatics analysis of sequence data
- Data interpretation

# LIBRARY PREPARATION

Library = fragments of DNA that have been prepared for amplification and sequencing


Genomic DNA fragment 150-600bp in length

The size of the fragment is called the “**insert size**”








# SEQUENCE ENRICHMENT OPTIONS



# PCR VS. SEQUENCE CAPTURE




# SEQUENCE CAPTURE



- Up to 80% enrichment for the targeted DNA
- 120 bp “baits” bind to DNA and magnetic beads
- Unbound DNA is discarded
- Baits are digested



# METHODOLOGIES FOR PERFORMING PCR BASED ENRICHMENT

- Uniplex PCR: 1 reaction = 1 amplicon
- Multiplex PCR: 1 reaction = 10-50 amplicons
- Droplet PCR: 1 reaction = 4,000 amplicons
- Microfluidics PCR
  - Multiplex 10 PCR/well
  - Can simultaneously amplify 480 amplicons



# CLUSTER GENERATION AND SEQUENCING


- Oligonucleotides attached to flow cell hybridize to the adaptors
- Individual DNA library fragments are immobilized






# CLUSTER GENERATION AND SEQUENCING

- Starting DNA template concentration is crucial to avoid overcrowding of clusters
- Each unique DNA molecule undergoes “bridge amplification”




# CLUSTER GENERATION AND SEQUENCING

- Simultaneous generation of millions of clusters (“colonies”)
- One cluster:
  - Derives from a single parent DNA molecule
  - Made up of ~1000 identical copies
  - Unique
  - Physically isolated from other clusters





# SEQUENCING BY SYNTHESIS

- All clusters are sequenced in parallel, one base at a time
- Fluorescently tagged nucleotides compete for next space
- Fluorescent tag blocks addition of more than 1 nucleotide per round
- Each round
  - Addition of one base
  - Laser excitation -> fluorescence
  - One “base” read from each cluster
  - Removal of fluorescent tag






# SEQUENCING BY SYNTHESIS



# BIOINFORMATICS ANALYSIS


## Phase 1: nGS data processing

Typically by lane




## Phase 2: variant discovery and genotyping

Typically multiple samples simultaneously but can be single sample alone



## Phase 3: integrative analysis



# INHERITED DISEASES VS. ONCOLOGY

	INHERITED DISEASES	ONCOLOGY
NUMBER OF GENES	Range: 1-150	Range: 1-50
MUTATION DISTRIBUTION	Across entire gene	Hotspot mutations/deletions /translocations
MINIMUM COVERAGE	20X	~250X – 1000X
STRUCTURAL VARIATION	Existing methods work well	Sensitivity depends on tumor percentage
TISSUE TYPES	Blood	Blood, Fresh frozen, FFPE, cytology
INPUT DNA	3 µg	5 ng – 1µg
TURNAROUND TIMES	4-6 weeks	5-7 days
COST	Range: \$1000 - \$10,000	Range: \$400 - \$1000

# TESTING SCHEME FOR INHERITED DISEASES

**SAMPLE  
PREPARATION**

**TARGET CAPTURE**



**SMALL PANELS**

RAPID TURNAROUND TIMES  
(5-7 DAYS)




**LARGE PANELS**

EXOME SEQUENCING  
SLOW TURNAROUND TIMES  
(8-10 WEEKS)




**SEQUENCING**

# INHERITED DISEASES:NGS TESTING



# NGS EXPERIENCE AT MDL

- Have offered NGS testing for over 130 Mendelian disorders since August 2012
  - We have tested 300 samples
  - We have detected mutations in approximately 30% of all samples tested
  - Mutation detection rate is dependent on clinical diagnosis
    - Mutation identified in 80% of inherited thrombophilias
    - Mutation identified in 25% of ataxias
    - No mutations identified in disorders of sexual development



# GENOMIC REGIONS THAT ARE PROBLEMATIC FOR NGS

**TRI-ALLELIC SNPS: 1  
DISCREPANCY (2%)**




**PSEUDOGENES/ HOMOLOGOUS  
REGIONS: 16 DISCREPANCIES (31%)**

**52 UNIQUE DISCREPANCIES IN CODING REGION OF DNA**

**GC RICH REGIONS :  
18 DISCREPANCIES  
(35%)**


**POLYMORPHIC  
REPEAT REGIONS:  
7 DISCREPANCIES  
(13%)**

**MISMAPPED LARGE  
INDELS :  
10 DISCREPANCIES  
(19%)**



# TESTING SCHEME FOR ONCOLOGY

**DNA/RNA EXTRACTED FROM SUBMITTED TISSUE**






# ONCOLOGY:NGS TESTING

- **Lung Cancer Panel**
- Somatic mutation testing
  - *KRAS (NRAS/HRAS)*
  - *EGFR*
  - *BRAF*
  - *PIK3CA*
  - *ERBB2*
  - *MET*
  - *TP53*
  - *AKT1*
  - *MAP2K1*
  - *EGFRvIII* (RT-PCR assay)
- Translocation
  - *ALK* (EML4-ALK, but other partners up to 20)
  - *ROS* (up to 7 partners)
  - *KIF5B/RET*
  - *CCDC6/RET* (aka RET/PTC1)
- Amplification
  - *EGFR*
    - *MET*
    - *MAPK1* (p42/ERK2)
    - *FGFR1*
    - *FGFR2*
- **Gastrointestinal Cancer Panel**
- Somatic mutation testing
  - *EGFR*
  - *KRAS (HRAS/NRAS)*
  - *BRAF*
  - *PIK3CA*
  - *TP53*
  - *ERBB2*
  - *MET*
  - *KIT*
  - *PDGFRA*
  - *AKT1*
  - *PTEN*
  - *APC*
- Amplification/Deletion
  - *ERBB2*
  - *IGF2* (11p15.5)
  - *PTEN*
  - *MDM2*
  - *EGFR* (rare)


# DETECTION OF EGFR MUTATIONS

- Targeted sequencing of exons 18-21 (visualizing 18-19)
- Input DNA: 50 ng of DNA for two lung cancer specimens
- Specimen 1 (exon 19 deletion)
  - 80% tumor
- Specimen 2 (L747P mutation due to sequential T>C mutations)
  - 70% tumor


# EGFR EXON 19 DELETION



# EGFR L747P MUTATION



# BRAF V600E MUTATION: MINIMUM DNA INPUT



# TECHNICAL ISSUES WITH NGS IMPLEMENTATION

- Bioinformatics methods for sequence alignment keep undergoing rapid improvements
  - Need to update bioinformatics pipeline at frequent intervals
- Structural genetic variation:
  - Optimal algorithms for detection of copy number variation remain unclear
- Several regions with inadequate coverage
  - Backup Sanger sequencing/alternative methodology necessary for several exons in the context of inherited disorders
  - Sensitivity to detect somatic mutations will not be the same in all the analyzed regions

# CLINICAL ISSUES WITH NGS IMPLEMENTATION


- Interpretation of clinical significance of many variants is unclear
  - Communication of these results to the clinician is problematic
  - Often results in additional testing of family members to determine clinical significance of a particular variant
- Incidental genetic findings need to be reported and appropriate clinical follow up procedures need to be in place

# OTHER ISSUES WITH NGS IMPLEMENTATION

- High upfront costs for test validation
  - Substantial reagent costs
- High sequencing run costs
  - Need to batch samples to reduce assay costs
    - Need to offer a large test menu to increase sample volume
  - Limited ability to repeat samples
    - Robustness of assays need to be adequately validated




# FUTURE APPLICATIONS OF NGS



# ACKNOWLEDGEMENTS

- **FUNDING SOURCES**
  - Institute for Translational Neuroscience
  - Biomedical Genomics Center
- **BIOMEDICAL GENOMICS CENTER**
  - Kenneth Beckman, Karina Bunjer, Adam Hauge, Archana Deshpande
- **BIOINFORMATICS CORE FACILITY**
  - Kevin Silverstein, Getiria Onsongo, Jesse Erdmann
- **MOLECULAR DIAGNOSTICS LABORATORY**
  - Matt Bower, Teresa Kemmer, Matt Schomaker, Sophia Yohe, Jon Wilson, Michael Spears, Andrew Nelson
- **FAIRVIEW**
  - Klint Kjeldahl, Karin Libby


THE MORE YOU KNOW, THE HARDER IT IS TO TAKE DECISIVE ACTION.



9/21

WATSON

ONCE YOU BECOME INFORMED, YOU START SEEING COMPLEXITIES AND SHADES OF GRAY.



© 1993 Waterson/Dell, by Universal Press Syndicate

YOU REALIZE THAT NOTHING IS AS CLEAR AND SIMPLE AS IT FIRST APPEARS. ULTIMATELY, KNOWLEDGE IS PARALYZING.



BEING A MAN OF ACTION, I CAN'T AFFORD TO TAKE THAT RISK.

YOU'RE IGNORANT, BUT AT LEAST YOU ACT ON IT.

